

<b>IHE Integration Statement</b>		<b>Date</b>	November 30, 2010
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Three Palm Software	WorkstationOne	1.5	
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Mammography Image Profile	Image Display	None	
	Evidence Creator	None	
	Print Composer	None	
Reporting Workflow Profile	Report Creator	None	
Key Image Note Profile	Image Display	None	
	Evidence Creator	None	
Mammography (Acquisition) Workflow Profile	Image Display	None	
	Evidence Creator	None	
Scheduled Workflow Profile	Image Display	None	
	Evidence Creator	None	
Internet address for TPS's IHE information: <a href="http://threepalmsoft.com/Support">threepalmsoft.com/Support</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
DICOM	<a href="#">WorkstationOneDICOMConformanceStatement.pdf</a>		
<b>Links to General Information on IHE</b>			
North America	<a href="http://www.ihe.net">www.ihe.net</a> <a href="http://wiki.ihe.net">wiki.ihe.net</a>		
Europe	<a href="http://www.ihe-europe.org">www.ihe-europe.org</a>		