

IHE Integration Statement		Date	January 2, 2015
Vendor	Product Name	Version	
Three Palm Software	WorkstationOne	1.7.2	
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image Profile	Image Display	None	
	Evidence Creator	None	
	Print Composer	None	
Digital Breast Tomosynthesis (DBT) (Trial status)	Image Display	Key Images, User Annotation	
	Evidence Creator	None	
	Print Composer	None	
Reporting Workflow Profile	Report Creator	None	
Key Image Note Profile	Image Display	None	
	Evidence Creator	None	
Mammography (Acquisition) Workflow Profile (Trial status)	Image Display	None	
	Evidence Creator	None	
Scheduled Workflow Profile	Image Display	None	
	Evidence Creator	None	
Internet address for TPS's IHE information: threepalmsoft.com/support			
Links to Standards Conformance Statements for the Implementation			
DICOM	WorkstationOneDICOMConformanceStatement.pdf		
Links to General Information on IHE			
North America	www.ihe.net wiki.ihe.net		
Europe	www.ihe-europe.org		